

Inspect unit, diagnose and isolate problem

11/18/2012 11:29 AM - tin

Status: Closed	Start date: 11/18/2012
Priority: Normal	Due date:
Assignee: tin	% Done: 100%
Category:	Estimated time: 0.00 hour
Target version:	
Description	

History

#1 - 11/18/2012 02:24 PM - tin

- % Done changed from 0 to 80

Bad contact somewhere, if assemble with case - UI go crazy.
Open case: Functions ok, self-test PASS, outputs ok.
Need cleaning board and replace DC Fan.

#2 - 06/03/2013 11:27 PM - tin

- Status changed from In Progress to Resolved

- % Done changed from 80 to 100

Unit fixed by replacing ICs from donor board.
Enabled to use in homelab.

#3 - 01/04/2014 05:11 PM - tin

- Status changed from Resolved to Closed